


RESTORATIVE



MAKE PEOPLE SMILE


This simple yet powerful idea drives every decision we make at Asa Dental. Whenever we choose materials, technologies, designs and manufacturing processes, we make sure they will help dental professionals make people smile, anywhere in the world. Our restorative instruments are the result of precise, careful planning and design which is rendered concrete through advanced technologies and skilled hands. This is how we manufacture durable excavators, precise carvers, effective burnishers, innovative isolation instruments and all the other products that allow practitioners to deliver aesthetically and functionally superior results.





TOOTH NUMBERING

Maxillary teeth	Patterns
18 17 16	4 W4 8 W8 8A W8A 12A 14 14T W14 14A W14A 26 51 205 202
15 14	0 W0 1 1T W1 2A 2AT W2A 29
13 12 11 21 22 23	9 9T 212 210
24 25	0 W0 1 1T W1 2A 2AT W2A 29
26 27 28	4 W4 8 W8 8A W8A 13A 14 14T W14 14A W14A 26 51
Mandibular teeth	Patterns
38 37 36	3 W3 7 W7 8A W8A 13A 14 14T W14 14A W14A
35 34	0 W0 2A 2AT W2A 2 2T W2 29
33 32 31 41 42 43	00 9 9T 212 210
44 45	0 W0 2 2T W2 2A 2AT W2A 29
46 47 48	3 W3 7 W7 12A 14 14T W14 14A W14A

CLASSIFICATION




DENTAL DAM CLAMPS Anterior



Premolars

for small, elongated or narrow teeth



for general use, with adjustable wings that adapt to the gingival curvature



for premolars, with large necks



Premolars

for lower premolars, general use



Molars

for square lower molars



for molar roots



for round, partially broken, undersized or slanted molars

Molars

for small lower molars



for upper molars



Molars

for small molars with round collars



4
winged
3052-4



W4
wingless
3052-W4



for round upper molars



8
winged
3052-8



W8
wingless
3052-W8



for partially broken molars



12A
winged, serrated
3052-12A



13A
winged, serrated
3052-13A



for large, partially broken, irregularly shaped or slanted molars



14A
winged
3052-14A



W14A
wingless
3052-W14A



for large upper molars



202
winged
3051-202



205
winged
3051-205



DENTAL DAM INSTRUMENTS

AsaLady

products ergonomically designed for women's hands:

- smaller and lighter handles
 - smaller finger holes
 - reduced size
- code starts with "W" or "SW"



latex dam
vacuum packed
36 pz / confezione



15x15 cm (6"x 6")
3045

Punch forceps



Ainsworth
17 cm
3032



Ainsworth
17 cm
W3036



Ivory
16.2 cm
3033

DENTAL DAM INSTRUMENTS

Rubberdamframe
stainless steel




10x10 cm
3040-1

Dental dam placed on an upper molar



Rubberdamframe
sterilizable plastic,
x-ray compatible



3041-1

Clamp forceps



Ivory lightweight
17 cm
3030



Ivory lightweight
16 cm
W3035



Brewer
17.5 cm
3031



Stoke
17.5 cm
3034

INSTRUMENT KIT FOR DENTAL DAMS




S3000

Dental dam instrument kit manufactured according to the instructions of Prof. Guido Goracci, Tenured Professor at the Department of Conservative Dentistry at La Sapienza University in Rome.

3040-1	stainless steel dental dam frame	Brewer clamp	
3031	forceps	Ivory punch	
3033	forceps	wingless clamp	
3050-60	for premolars	wingless clamp	60
3051-26	clamp for molars	wingless	26
3051-51	clamp for premolars	winged	51
3051-29	winged clamp for molars	winged	29
3051-202	winged posterior	clamp	202
3051-205	winged anterior	clamp	205
3051-212	stainless steel dental dam		212
3051-210	clamp tray		210
2731-1	stainless steel tray		
	stainless steel lid		
2720-L	stainless steel	insert	
2720-U	frame		
2730-3			

Insert frame for dental dam tools



2730-3


INSTRUMENT KIT FOR DENTAL DAMS



S2731-1

stainless steel dental dam clamp tray with 8 clamps

3050-60	wingless clamp for	60
3051-26	premolars wingless clamp	26
3051-51	for molars wingless clamp for	51
3051-29	molars wingless clamp for	29
3051-202	premolars wingless clamp for molars	202
3051-205	winged clamp for molars	205
3051-212	winged posterior clamp	212
3051-210	winged anterior clamp	210
2731-1	stainless steel dental dam clamp tray	



S2731-12

stainless steel dental dam clamp tray with 12 clamps

3051-212	winged posterior clamp	212
3052-00	winged premolar clamp	00
3052-0	winged premolar clamp	0
3052-1	winged premolar clamp	1
3052-2A	winged premolar clamp	2A
3052-3	winged clamp for molars	3
3052-7	winged clamp for molars	7
3052-8	winged clamp for molars	8
3052-8A	winged molar root clamp	8A
3052-12A	serrated and winged molar clamp	12A
3052-13A	serrated and winged molar clamp	13A
3052-14	winged clamp for molars	14
2731-1	stainless steel dental dam clamp tray	



2731-1

stainless steel dental dam clamp tray (without clamps)


HANDLE

SINGLE-ENDED

Asa35
(stainless steel)
ø 3.5 mm

Asa60
(stainless steel)
ø 6.0 mm

Octagonal
(stainless steel)
ø 5.5 mm



DOUBLE-ENDED

Asa60
(stainless steel)
ø 6.0 mm

Asa Light MV
(stainless steel, hollow)
ø 8.0 mm

Leonardo ML
(stainless steel)
ø 7.0 mm

Octagonal
(stainless steel)
ø 5.5 mm

Alu80
(aluminum, hollow)
ø 8.0 mm

AsaWave RC
(plastic resin in PPSU with a steel core)
ø 7.5 / 14.5 mm



HOW TO USE THE INFORMATION


- 1/2 Wedelstaedt → model + scientific name
- MV1604-02 → Asa Dental code
- Asa Light MV → handle option

AsaLady

products ergonomically designed for women's hands:

- smaller and lighter handles
- smaller finger holes
- reduced size

code starts with "W" or "SW"



EXCAVATORS
Curettes




Set of 5



Model	Set Code	Handle	Asa
0 · 1 · 2 · 3 · 4	SMV1700	Light	MV
0 · 1 · 2 · 3 · 4	S1700	Asa60	
4	S1700-1	Asa60	
0	S1700-2	Asa60	
1	S1700-3	Asa60	
2	S1700-4	Asa60	
3			
4			

EXCAVATORS
Curettes



EXCAVATORS
Oval currettes



1.2 mm 1.5 mm 2 mm



17W
1704-17W
Asa60
MV1704-17W
Asa Light MV

18W
1705-18W
Asa60
MV1705-18W
Asa Light MV

19W
1706-19W
Asa60
MV1706-19W
Asa Light MV

EXCAVATORS
ASH



1 mm 1.4 mm



243
1711-243
Asa60
ML1711-243
Leonardo ML

G2
1711-G2
Asa60
ML1711-G2
Leonardo ML

EXCAVATORS
SS. white type



Set of 6

Model 5 · 14 · 17 ·
18 · 19 · 20 5 · 14 ·
17 · 18 · 19 · 20

Set Code	Handle
S1710-6	Asa60
SML1710-6	Leonardo ML

1 mm 1.2 mm 1 mm 2 mm



5
1710-5
Asa60
ML1710-5
Leonardo ML

14
1710-14
Asa60
ML1710-14
Leonardo ML

17
1710-17
Asa60
ML1710-17
Leonardo ML

18
1710-18
Asa60
ML1710-18
Leonardo ML

19
1710-19
Asa60
ML1710-19
Leonardo ML

20
1710-20
Asa60
ML1710-20
Leonardo ML

EXCAVATORS
Darby-Perry

0.8 mm 0.8 mm 1 mm 1 mm 1 mm 2 mm

2 Darby-Perry
1711-2
Asa60
ML1711-2
Leonardo ML

11 Darby-Perry
1711-11
Asa60
ML1711-11
Leonardo ML

12 Darby-Perry
17 11 -12
Asa60
ML1711-12
Leonardo ML

21 Darby-Perry
1711-21
Asa60
ML1711-21
Leonardo ML

1.2 mm 1.6 mm 1.6 mm 2 mm

1 Darby-Perry
ML1711-1
Leonardo ML
17 11 -1
Asa60L

22 Darby-Perry
1711-22
Asa60
ML1711-22
Leonardo ML

31 Darby-Perry
1711-31
Asa60
ML1711-31
Leonardo ML

32 Darby-Perry
1711-32
Asa60
ML1711-32
Leonardo ML

41 Darby-Perry
1711-41
Asa60
ML1711-41
Leonardo ML

42 Darby-Perry
ML1711-42
Leonardo ML
1711-42
Asa60

CAVITY PREPARATION INSTRUMENTS
Gingival margin trimmers

2 Black/27
mesial
1601-2
Asa60
MV1601-2
Asa Light MV

3 Black/28
distal
1601-3
Asa60
MV1601-3
Asa Light MV

1 Black/26
distal
1601-1
Asa60
MV1601-1
Asa Light MV

4 Black/29
mesial
1601-4
Asa60
MV1601-4
Asa Light MV

1 Black/77-78
distal
1 6 02-1
Asa60
MV1602-1
Asa Light MV

2 Black/79-80
mesial
1 6 02-2
Asa60
MV1602-2
Asa Light MV

CAVITY PREPARATION INSTRUMENTS
Hatchets



- | | |
|--|--|
| 1 Black/53-54
1603-1
Asa60
MV1603-1
Asa Light MV | 2 Black/51-52
1603-2
Asa60
MV1603-2
Asa Light MV |
|--|--|

CAVITY PREPARATION INSTRUMENTS
Chisels



- | | | |
|--|--|--|
| 1/2 Wedelstaedt
1604-02
Asa60
MV1604-02
Asa Light MV | 3/4 Wedelstaedt
1604-04
Asa60
MV1604-04
Asa Light MV | 5/6 Wedelstaedt
1604-06
Asa60
MV1604-06
Asa Light MV |
|--|--|--|

AMALGAM CARRIERS



AMALGAM WELLS
cold sterilization only



Ø 3.5 cm
h. 3.5 cm
1030-1



Ø 2 cm
h. 4 cm
1030-2



Ø 2 cm
h. 2.5 cm
1030-3



Ø 1.5 mm
Ø 2 mm
medium/small
20 cm
1015-1



Ø 2.5 mm
Ø 2 mm
medium/large
20 cm
1015-2



Ø 3 mm
Ø 2 mm
medium/jumbo
20 cm
1015-3



Ø 1.5 mm
small
15.5 cm
1011-1



Ø 2 mm
medium
15.5 cm
1011-2



Ø 2.5 mm
large
15.5 cm
1011-3



Ø 3 mm
jumbo
15.5 cm
1011-4

AMALGAM CARRIERS




CEMENT SPATULAS

Set of 5



Model	Set Code
1 2 3	S0903-1
	S0903-2
	S0903-3




HEIDEMANN SEPARATING SPATULAS



HEIDEMANN SEPARATING SPATULAS

Set of 5

Model	Set Code	Handle
1 2 3	S 11 02-	Asa60
	1	Asa60
	S1102-2	Asa60
	S1102-3	



HEIDEMANN SEPARATING SPATULAS with titanium nitride-coated tips



AMALGAM FILLING INSTRUMENTS
with titanium nitride-coated tips

Set of 5



Model	Code	Set Code	Handle
1 2 1	1200T-	S1200T	Asa60
1 2	1	S1200T	Asa60
	1202T-	S1200T	Asa60
	2 1 2 0	S1200T	Asa60
	4T-1	S1200T	Asa60



1200T-2




1.2 mm
1
1200T-1
Asa60

AMALGAM FILLING INSTRUMENTS

Set of 5



Model	Set Code	Handle
1 2 3	S1200-1	Asa60
4 5	S1200-2	Asa60
	S1200-3	Asa60
	S1200-4	Asa60
	S1200-5	Asa60



AMALGAM FILLING INSTRUMENTS

1.2 mm serrated
1.8 mm serrated
2.3 mm serrated
3 mm serrated
1.2 mm non-serrated
1.8 mm non-serrated

2 mm
2.6 mm
3.2 mm
1.4 mm
2 mm

2
1202-2
Asa60
MV1202-2
Asa Light MV
RC1202-2
AsaWave RC

3
1202-3
Asa60
MV1202-3
Asa Light MV

4
1202-4
Asa60
MV1202-4
Asa Light MV

1L
1202-1L
Asa60
MV1202-1L
Asa Light MV

2L
1202-2L
Asa60
MV1202-2L
Asa Light MV

2.3 mm non-serrated
3 mm non-serrated
1 mm
1.4 mm
1.8 mm

1.4 mm
2.6 mm
3.2 mm
1.2 mm
1.6 mm
2 mm

1
1202-1
Asa60
MV1202-1
Asa Light MV

3L
1202-3L
Asa60
MV1202-3L
Asa Light MV

4L
1202-4L
Asa60
MV1202-4L
Asa Light MV

1
1203-1
Asa60
MV1203-1
Asa Light MV

2
1203-2
Asa60
MV1203-2
Asa Light MV

3
1203-3
Asa60
MV1203-3
Asa Light MV

AMALGAM FILLING INSTRUMENTS

2.4 mm
1.4 mm
2.4 mm
1.2 mm
1.8 mm

2.8 mm
0.8 mm
1.4 mm
0.8 mm
1.5 mm

2 Ward
1205-2
Asa60
MV1205-2
Asa Light MV

0T Tanner
1211-0T
Asa60
MV1211-0T
Asa Light MV

1T Tanner
1211-1T
Asa60
MV1211-1T
Asa Light MV

2T Tanner
1211-2T
Asa60
MV1211-2T
Asa Light MV

3T Tanner
1211-3T
Asa60
MV1211-3T
Asa Light MV

1 mm serrated
1.5 mm
3 mm
2.6 mm
1.8 mm

1.4 mm
1 mm
2 mm
1 mm
2.6 mm

01 Marquette
1212-01
Asa60
MV1212-01
Asa Light MV

H1 Hollenback
1213-H1
Asa60
MV1213-H1
Asa Light MV

H2 Hollenback
1213-H2
Asa60
MV1213-H2
Asa Light MV

H3 Hollenback
1213-H3
Asa60
MV1213-H3
Asa Light MV

H4 Hollenback
1213-H4
Asa60
MV1213-H4
Asa Light MV

2 mm
1 Ward
1205-1
Asa60
MV1205-1
Asa Light MV

AMALGAM FILLING INSTRUMENTS
Markley





Set of 6









Model 1 · 2 · 3
· 4 · 5 · 6 1 · 2 ·
3 · 4 · 5 · 6

Set Code
S1400
SMV1400

Handle
Asa60
Asa Light MV

<p>0.4 mm</p>  <p>0.6 mm</p> <p>1 Markley 1400-1 Asa60 MV1400-1 Asa Light MV</p>	<p>0.5 mm</p>  <p>1 mm</p> <p>2 Markley 1400-2 Asa60 MV1400-2 Asa Light MV</p> <p>1.6 mm</p> <p>4 Markley 1400-4 Asa60 MV1400-4 Asa Light MV</p>	<p>0.4 mm</p>  <p>1.2 mm</p> <p>3 Markley 1400-3 Asa60 MV1400-3 Asa Light MV</p> <p>0.5 mm</p> <p>5 Markley 1400-5 Asa60 MV1400-5 Asa Light MV</p>	<p>3 mm</p>  <p>1.6 mm</p> <p>6 Markley 1400-6 Asa60 MV1400-6 Asa Light MV</p>
---	--	--	---

AMALGAM FILLING INSTRUMENTS
Foc

 <p>AsaLady</p> <p>0 Foc ML1401-0/WML1401-0 Leonardo ML 1401-0 Asa60</p>	<p>serrated</p>  <p>0S Foc 1401-0S Asa60 ML1401-0S Leonardo ML</p> <p>serrated</p>  <p>1S Foc 1401-1S Asa60 ML1401-1S Leonardo ML</p>	<p>serrated</p>  <p>1 Foc 14 01 -1 Asa60 ML1401-1 Leonardo ML</p> <p>AsaLady</p>  <p>2 Foc 1401-2 Asa60 ML1401-2 Leonardo ML</p>	 <p>2S Foc 1401-2S Asa60 ML1401-2S Leonardo ML</p>
--	---	---	---

FILLING INSTRUMENTS
Burnishers

2.2 mm
2.8 mm
1.2 mm
1 mm
1.4 mm
2 mm
3.2 mm
2 Westcott
1204-2
Asa60
MV1204-2
Asa Light MV
2 mm
2 Ladmore
1214-2
Asa60
MV1214-2
Asa Light MV
1.8 mm
3 Ladmore
1214-3
Asa60
MV1214-3
Asa Light MV
2.5 mm
18 Ball
1215-18
Asa60
MV1215-18
Asa Light MV
RC1215-18
● AsaWave RC
1 mm
26-27S Ball
1215-27S
Asa60
MV1215-27S
Asa Light MV
RC1215-27S
● AsaWave RC
1.6 mm
2.5 mm
3 mm
1.3 mm
2.5 mm
1 Westcott
1204-1
Asa60
MV1204-1
Asa Light MV
3 mm
27-29 Ball
1215-29
Asa60
MV1215-29
Asa Light MV
1.3 mm
A
1214-A
Asa60
3.4 mm
W3
1215-W3
Asa60
NY U/B2
1220-B2
Asa60
MV1220-B2
Asa Light MV
2/31
1221-31
Asa60
MV1221-31
Asa Light MV

PLASTIC FILLING INSTRUMENTS
with titanium nitride-coated tips

1 mm
1.6 mm
2 mm
2 mm
2.9 mm
1.2 mm
2 mm
2.6 mm
0.1 mm
0.4 mm
1
1208T-1
Asa60
2
1208T-2
Asa60
3
1208T-3
Asa60
2
1209T-2
Asa60
2
1210T-2
Asa60

PLASTIC FILLING INSTRUMENTS

1.2 mm, 1.6 mm, 1.2 mm, 1.6 mm, 1 mm, 2.5 mm

0.1 mm, 2.2 mm, 2.7 mm, 2 mm

0.7 mm, 2.5 mm, 1 mm, 2 mm

0.1 mm, 2.5 mm, 3 mm, 0.1 mm, 0.1 mm, 0.1 mm

1 1206-1 Asa60 MV1206-1 Asa Light MV

2 1206-2 Asa60 MV1206-2 Asa Light MV

1 1207-1 Asa60 MV1207-1 Asa Light MV

2 1207-2 Asa60 MV1207-2 Asa Light MV

1 1208-1 Asa60 MV1208-1 Asa Light MV

1 1209-1 Asa60 MV1209-1 Asa Light MV

2 1209-2 Asa60 MV1209-2 Asa Light MV

3 1209-3 Asa60 MV1209-3 Asa Light MV

PLASTIC FILLING INSTRUMENTS

2 1210-2 Asa60 MV1210-2 Asa Light MV RC1210-2 AsaWave RC

3 1210-3 Asa60 MV1210-3 Asa Light MV RC1210-3 AsaWave RC

1 Woodson 1216-1 Asa60 MV1216-1 Asa Light MV

2 Woodson 1216-2 Asa60 MV1216-2 Asa Light MV

1 1210-1 Asa60 MV1210-1 Asa Light MV RC1210-1 AsaWave RC

3 Woodson 1216-3 Asa60 MV1216-3 Asa Light MV

11 1217-1 1 Asa60 MV1217-1 1 Asa Light MV

6 1218-6 Asa60 MV1218-6 Asa Light MV

179 1219-179 Asa60 MV1219-179 Asa Light MV

MODELING INSTRUMENTS

1 Ward
1501-1
Asa60
MV1501-1
Asa Light MV

1S Ward
1501-1S
Asa60
MV1501-1S
Asa Light MV

1/2 Hollenback
1500-02
Asa60
MV1500-02
Asa Light MV

2 Ward
1501-2
Asa60
MV1501-2
Asa Light MV

3 Hollenback
1500-3
Asa60
MV1500-3
Asa Light MV

3S Hollenback
1500-3S
Asa60
MV1500-3S
Asa Light MV

1 Frahm
15 02-1
Asa60
M V 15 02-1
Asa Light MV

MODELING INSTRUMENTS

2 USC
1503-2
Asa60

1 Interproximal
1504-1
Asa60

2 Interproximal off-angle
1504-2
Asa60

SHO -A
1505-A
Asa60

3 -SHO
1505-3
Asa60

104
1506-104
Asa60

8 Wiland
1507-8
Asa60
MV1507-8
Asa Light MV

MODELING INSTRUMENTS

Instruments for retraction cords



- 1 cleoid- discoïd 1600-1 Asa60 MV1600-1 Asa Light MV
- 2 cleoid- discoïd 1600-2 Asa60 MV1600-2 Asa Light MV
- 3 cleoid- discoïd 1600-3 Asa60 MV1600-3 Asa Light MV
- 4 cleoid- discoïd 1600-4 Asa60 MV1600-4 Asa Light MV

- 1 1310-1 Asa60
- 2 1310-2 Asa60

OPERATIVE INSTRUMENTS

Instrument for the application of tinctures and medicaments

Placement instruments



Gauze packers




- 1300 Asa35
- 1301 Asa60 MV1301 Asa Light MV
- 1302 Asa60 MV1302 Asa Light MV


- 1 Luniatschek 2 mm 1355-1 Octagonal
- 2 Luniatschek 3 mm 1355-2 Octagonal

CROWN INSTRUMENTS

Crown remover
complete with double-ended tips



Crown remover
complete with single-ended tips



CROWN INSTRUMENTS

Crown removers and spreaders



CROWN INSTRUMENTS
Decimal calipers



measuring range: 0-10 mm
(indicative measures)



Iwanson
for wax
2550-1

measuring range: 0-100 mm
(indicative measures)



Boley
14
cm
2555

measuring range: 0-10 mm
(indicative measures)



Iwanson
for metal
2550-2

Articulating paper forceps




Miller
15
cm
0670




15 cm
0671

CROWN INSTRUMENTS
Telescopic crown pliers



A	4.4 mm
B	5.0 mm
C	3.4 mm
D	3.8 mm



A	3.3 mm
B	4.0 mm
C	3.8 mm
D	3.0 mm



Böhm
14.5 cm
2500-2



Böhm
11.5 cm
2500-1

Copper ring removing pliers



Furrer
straight
16 cm
2331-1



Furrer
curved
16 cm
2331-2

HEMOSTATS



1 Halstead-Mosquito
straight
12 cm
0500-1



1 Kelly
straight
14 cm
0510-1



1 Micro-Mosquito
straight
10 cm
0501-1

SCISSORS
for gold leaf



1
straight
12 cm
0300-1

UNIVERSAL SCISSORS
with tungsten carbide blades



TC Universal
12 cm
5700TC

UNIVERSAL SCISSORS
one serrated blade



Universal
12 cm
5700

CROWN SCISSORS



1
straight
one serrated blade
10.5 cm
5600-1



2
curved
one serrated blade
10.5 cm
5600-2



1
straight
10.5 cm
5601-1



2
curved
10.5 cm
5601-2

MATRIX RETAINERS



MATRICES strip

for matrix retainers

- 3070-5
- 3071
- 3072
- 3073-1
- 3073-2
- 3073-3

for matrix retainers

- 3070-6
- 3070-7
- 3071
- 3072
- 3073-1
- 3073-2
- 3073-3

for matrix retainers

- 3070-7
- 3071
- 3072
- 3073-1
- 3073-2
- 3073-3



5 mm x 100 cm 3082-5



6 mm x 100 cm 3082-6



7 mm x 100 cm 3082-7

Tofflemire 12 pcs / pack

12 pcs / pack assorted:



1 Tofflemire 3081-1



2 Tofflemire 3081-2



3 Tofflemire 3081-3

- 3081-1
- 3081-2
- 3081-3



Tofflemire 3081

COMPOSITE INSTRUMENTS
Molars and Pre-molars

Molars, Pre-molars, Incisors and Canines



Asa Com 1/1B
two tips in the same shape but different sizes, for filling and modeling occlusal anatomy in premolars and molars
Universal

Asa Com 2/2B
same tip but with two different angles to model occlusal ridges in premolars and molars
Universal

Asa Com 3/3B
highly versatile condenser with small and large semi-spherical tip, for placement and contouring valid for: classes 1, 2, 3, 4 and 5

Asa Com 4/4B
highly versatile filling instrument with different size semi-spherical tips, for placement and contouring; can be used for: classes 1, 2, 3, 4 and 5

Asa Com 5/5B
versatile instrument for the application, layering and modeling of composite, with opposite (lateral and frontal) blades; can be used for: classes 1, 2, 3, 4 and 5

Asa Com 6/6B
filling and shaping instrument ideal for use in mesial and distal areas; can be used for: classes 2, 3 and 5

Asa Com 1/1B
5112-3
Asa60

Asa Com 2/2B
0704-16D2
Asa60

Asa Com 3/3B
1203-1
Asa60

Asa Com 4/4B
1203-3
Asa60

Asa Com 5/5B
1219-179
Asa60

Asa Com 6/6B
1504-2
Asa60

COMPOSITE INSTRUMENTS
Molars, Pre-molars, Incisors and Canines

Molars and Pre-molars



Asa Com 7/7B
versatile instrument ideal for gingival areas but also suitable for mesial and distal definition; can be used for: classes 1, 2, 3, 4 and 5

Asa Com 8/8B
instrument that combines a small tip for occlusal anatomy modeling and a narrow paddle for anterior areas
Universal

Asa Com 9/9B
condenser equipped with 2 different-size flat tips that can be used for: classes 1, 2, 3, 4 and 5

Com 10/10B
condenser equipped with 2 different-size flat tips that can be used for: classes 1, 2, 3, 4 and 5

Com 11/11B
small and large tip conical condenser
Universal

Asa Com 12/12B
universal condenser with opposite (lateral and frontal) tips
Universal

Asa Com 7/7B
5153-1
● Alu80

Asa Com 8/8B
1220-B2
Asa60

Asa Com 9/9B
1201-1L
Asa60

Asa Com 10/10B
1201-2L
Asa60

Asa Com 11/11B
14 01 -1
Asa60

Asa Com 12/12B
1213-H3
Asa60

COMPOSITE INSTRUMENTS

Molars, Pre-molars, Incisors and Canines

AsaWave composite instruments



Asa Com 14/14B
1210-1
Asa60
MV1210-1
Asa Light MV
RC1210-1
● AsaWave RC

Asa Com 13/13B instrument with opposite blades for composite contouring, removal and modeling
Universal

Asa Com 14/14B versatile instrument for composite application, layering and modeling, with opposite (lateral and frontal) tips
Universal

Asa Com 15/15B versatile instrument for composite application, layering and modeling, with opposite (lateral and frontal) tips; can be used for: classes 2, 3, 4 and 5

Asa Com 13/13B
5112-4
Asa60

Asa Com 15/15B
210-3
Asa60
MV1210-3
Asa Light MV
RC1210-3
● AsaWave RC



RC1510-2D
● AsaWave RC

RC1510-3D
● AsaWave RC

RC1510-1D
● AsaWave RC

RC1510-4D
● AsaWave RC

RC1510-5D
● AsaWave RC

COMPOSITE INSTRUMENTS



S1222-7



S1222-15



Asa Com 1/1B 5112-3
Asa Com 2/2B 0704-16D2
Asa Com 3/3B 1203 -1
Asa Com 4/4B 1203 -3
Asa Com 5/5B 1219-179
Asa Com 6/6B 1504-2
Asa Com 7/7B 5153 -1

Asa Com 1/1B 5112-3
Asa Com 2/2B 0704-16D2
Asa Com 3/3B 1203 -1
Asa Com 4/4B 1203 -3
Asa Com 5/5B 1219-179
Asa Com 6/6B 1504-2
Asa Com 7/7B 5153 -1
Asa Com 8/8B 1220 - B2
Asa Com 9/9B 1201 -1L
Asa Com 10/10B 1201-2L
Asa Com 11/11B 1401-1
Asa Com 12/12B 1213 - H3
Asa Com 13/13B 5112-4
Asa Com 14/14B 1210 -1
Asa Com 15/15B 1210 -3