



# Electrosurgical Instruments



PRODUCT CATALOGUE

# Contents



Reusable Electrodes	<b>3</b>
Single Use Bi Polar Forceps	<b>4</b>
Standard Reusable Bi Polar Forceps	<b>5</b>
Bayonet Reusable Bi Polar Forceps	<b>6</b>
Non-Stick Reusable Bi Polar Forceps	<b>7</b>
Bi Polar Cables	<b>8</b>
Standard Reusable Mono Polar Forceps	<b>9</b>
Mono Polar Cables	<b>10</b>
Finger Switching Pencils Reusable & Single Use	<b>11</b>
Foot Control Pencils Reusable	<b>12</b>
Reusable Adaptors	<b>13</b>
Screw Based Quivers SQUIVER™	<b>14</b>
Repair and Reinsulation Service	<b>15</b>

# Reusable Electrodes



## Blade Electrodes



**Standard Blade**

300 - 25 - 000



**Extended Blade**


300 - 25 - 075

## Ball Electrodes



**Standard Ball**

400 - 4 - 000



**Extended Ball**

400 - 4 - 075

## Needle Electrodes




**Micro Dissection (Tungsten)**

PMS - 1600



**Standard Needle**

500 - 25 - 000/.65



**Extended Needle**

500 - 25 - 075/.65

## Laparoscopic Hook Electrodes

4 mm Banana Plug



**Reusable**

802 - 330 - 5H

5 mm Ø

# Single Use Bi Polar Forceps

## Micro Jewellers

A210 - 140 - 05



## Adson (smooth) Adson (serrated)

A210 - 200 - 10    A210 - 200 - 10SR



## Cushing Straight

A210 - 230 - 10



## Jansen Bayonet


A210 - 370 - 10



Code	A210-	Description	Tip Size	Tip Style	Length
140-05	A210-	Micro Jewellers	0.5mm	Smooth	115mm
200-10	A210-	Adson Adson	1.0mm	Smooth	120mm
200-10SR		Cushing Straight	1.0mm	Serrated	120mm
A210-230-10		Jansen Bayonet	1.0mm	Smooth	150mm
A210-370-10			1.0mm	Smooth	200mm

## Forcep Tip Size Guide

Single Use Bi Polar Cable 2 mm  
Code: SW46211



# Standard Reusable Bi Polar Forceps



**McPherson Angled**

701 - 100A



**Micro Jewellers**

701 - 115S



**Adson**

701 - 150BIA



**Semken**

701 - 150BI



**Cushing**

701 - 175BI



**Potts-Smith**

701 - 200BI

Code	Description	Tip Size	Tip Style	Length
701-100A	McPherson Angled	0.5mm	Smooth	100mm
701-115S	Jewellers Straight	0.5mm	Smooth	115mm
701-115A	Jewellers Angled	0.5mm	Smooth	115mm
701-150BIA	Adson Adson	0.5 or 1.0mm	Smooth	120mm
701-150BIA/SR	Semken Cushing	0.5 or 1.0mm	Serrated	120mm
701-150BI	Potts-Smith	0.5 ,1.0 or 1.5mm	Smooth	150mm
701-175BI		0.5 ,1.0 or 1.5mm	Smooth	175mm
701-200BI		0.5 ,1.0 or 1.5mm	Smooth	200mm

# Bayonet Reusable Bi Polar Forceps



**Bayonet Straight**  
701 - 175BS



**Bayonet Angled**  
701 - 175BA



**Bayonet Straight**  
701 - 200BS



**Bayonet Angled**  
701 - 200BA



**Bayonet Irrigating**  
701 - 200BS/IR

Code	Description	Tip Size	Tip Style	Length
701-175BS	Bayonet Straight	0.5, 1.0, 1.5mm	Smooth	175mm
701-175BA	Bayonet Angled	0.5, 1.0, 1.5mm	Smooth	175mm
701-200BS	Bayonet Straight	0.5, 1.0, 1.5mm	Smooth	200mm
701-200BA	Bayonet Angled	0.5, 1.0, 1.5mm	Smooth	200mm
701-200BS/IR	Bayonet Irrigating	1.0mm	Smooth	200mm

# Non-Stick Reusable Bi Polar Forceps



**Micro Jewellers**  
701 - 115S / NS



**Adson**  
701 - 150BIA / NS



**Cushing Straight**  
701 - 150BI / NS



**Jansen Bayonet**  
701 - 200BS / NS

## Non-Stick Bipolar Forceps

- Significant reduction of tissue adhesion to tips
- Reduced tissue trauma
- Highly conductive tip greatly reduces heat

Our **Non-Stick** range of bipolar forceps offers significant advantages over standard forceps when used in coagulation procedures. The adhesion of “coagulated” tissue to the tips is substantially reduced due to the heat dispersion properties of the unique metal alloy tipped forceps; this alloy is the best heat conducting material available and the only material which is bio - compatible.

Code	Description	Tip Size	Tip Style	Length
701-115S/NS	Micro Jewellers	0.5mm	Smooth	115mm
701-150BIA/NS	Adson	1.0mm	Smooth	120mm
701-150BI/NS	Cushing Straight	1.0mm	Smooth	150mm
701-200BS/NS	Jansen Bayonet	1.0mm	Smooth	200mm

# Bi Polar Cables

## Reusable

### V-16

- Round instrument plug
- 3 metres
- Silicone construction
- To suit twin pin forceps



## Reusable

### SW46200A

- Angled instrument plug
- 3 metres
- Silicone construction
- To suit twin pin forceps



## Single use

### SW46211

- Angled instrument plug
- 3 metres
- To suit twin pin forceps
- Individually packed in Steri-peel pouches. (box of 100)



# Standard Reusable Mono Polar Forceps



**Adson**  
700 - 140AP



**Potts-Smith**  
700 - 150ST



**Potts-Smith**  
700 - 175ST



**Potts-Smith**  
700 - 200ST



**Bayonet Straight**  
700 - 200BS



**Bayonet Angled**  
700 - 200BA

**Note:** Available  
Terminals - 5/32  
screw - 2.4mm  
Female pin - 4mm  
Banana plug

Code	Description	Tip Size	Tip Style	Length
700-140AP	Adson Potts-	0.5, 1.0 mm	Serrated	120 mm
700-150ST	Smith Potts-Smith	0.5, 1.0, 1.5 mm	Serrated	150 mm
700-175ST	Potts-Smith	0.5, 1.0, 1.5 mm	Serrated	175 mm
700-200ST	Bayonet - Straight	0.5, 1.0, 1.5 mm	Serrated	200 mm
700-200BS	Bayonet - Angled	0.5, 1.0, 1.5 mm	Serrated	200 mm
700-200BA		0.5, 1.0, 1.5 mm	Serrated	200 mm

# Mono Polar Cables



## Reusable

**900-02MJ / 900 - 03P**

- 2.4mm pin to 8mm
- machine plug
- Reusable
- 3 metres



## Reusable

**900-02MJ / 900 - 055F**

- 4mm socket to 8mm
- machine plug
- Reusable
- 3 metres



## Reusable

**900-02MJ / 900 - 06ST**

- 5/32 screw fitting to
- 8mm machine plug
- Reusable
- 3 metres



## Single use

**SW47311**


- 4mm socket to 8mm
- machine plug
- Single use
- 3 metres
- Individually packed in steri-peel pouches (box of 100)

# Finger Switching Pencils

## Reusable & Single Use

**Reusable**

**PMS750**



**Standard Single Use**

**SW12200**



**Non-Stick Single use**

**SW12200-T**



Code	Description	Box Quantity
PMS750	Reusable finger switching pencil with 3 metre silicone cable	1
SW12200	Single use finger switching pencil with stainless steel blade electrode	100
SW12200-T	Single use finger switching pencil with non-stick blade electrode	100

# Foot Control Pencils

## Reusable



801 - 200 - S



801 - 150 - S

5/32" Screw fitting




801 - 200 - P



801 - 150 - P

Pin Fitting



801 - 200 - 4



801 - 150 - 4

4mm Banana Plug



PMS - 652



Foot-controlled pencil with 3.0 m cable and 8.0 mm plug.

Code 801-	Length	Length	5/32" Screw	Pin Fitting	4 mm
200-S	8"	200 mm	Yes	-	-
801-150-S	6"	150 mm	Yes	-	-
801-200-P	8"	200 mm	-	Yes	-
801-150P	6"	150 mm	-	Yes	-
801-200-4	8"	200 mm	-	-	Yes
801-150-4	6"	150 mm	-	-	Yes

# Reusable Adaptors

## Machine Active Plug Adaptor

To suit all 4 mm banana plugs.



900-01-LJA

## Machine Active Plug Adaptor

To take hand-switching pens to foot-switching



900-01-HS/FS

## Bi-Polar Machine Adaptors

To suit CONMED Generators



ALI-BPA-32

To suit ERBE Generators



ALI-BPA-22

# Screw Based Quivers

## SQUIVER™




### Screw off base enables easy cleaning

- Autoclaveable - standard steam
- Belled top for easy entry
- Screw base will not fall or drop off whilst in use
- Custom made for telescopes or extra long laparoscopic probes
- Spare caps and stainless steel rings are available

### Sizes available:

Code	Dimensions
DQS1.5	150mm
DQS2	200mm
DQS2.5	250mm
DQS3	300mm
DQS3.5	350mm
DQS4	400mm

# Repair and Reinsulation Service



**Before**



**After**

## It's the environmentally responsible choice.

Alimtype specialises in the repair, reinsulation and complete refurbishment of all types and brands of electro-surgical instruments.

We can restore your old, worn and damaged instruments to as new condition using only the highest quality materials and techniques, adhering to the most stringent standards.

Call us today.

**Narham Enterprises, a Sialkot-based Pakistani company, has been at the forefront of design, manufacture, and servicing of surgical instrumentation for over 8 years.**



## **Manufacturing / Surgical Supplies**

**Narham Enterprises manufactures and supplies a wide range of general and specialised surgical instruments, including:**

- General theatre instruments
- Laproscopy instruments
- Orthopaedic instruments
- Special technique needles, cass waldie etc
- Biopsy needles
- Custom made instruments
- Modifications and prototyping

## **Repairs / Servicing**

**Narham Enterprises offers a complete repair and maintenance service, which includes:**

- Sharpening - scissors, osteotomes, rongeurs
- Reblocking needle holders
- Reinsulation of diathermy instruments
- General repairs and maintenance of all instruments
- Super fast turn around time

### **How to order**

-  By Phone on **(03) 700048270**
-  By Fax to **(03) 700048270**
-  Email: **info@narhamentp.com**
-  Contact your local representative
-  Roras Road, Muzaffarpur, Sialkot  
51310 Pakistan.

**Call for competitive pricing**

 **(03) 700048270**