




CARTILAGE AND TENDON INSTRUMENTS

Cartilage and tendon instruments are specially designed for delicate handling and precise manipulation of soft tissues during orthopedic procedures. Our instruments ensure minimal trauma while providing accurate control for cutting, grasping, or repositioning cartilage and tendons. Crafted from high-quality surgical steel, they support surgeons in performing meticulous repairs with confidence and care.

COBB ELEVATORS

COBB SPINAL ELEVATORS, HEX HANDLE


All stainless steel and hexagonal handle. Sharp, 9-1/2" (24 cm).



CODE #	DESCRIPTION
NE49340	1/2" (13 mm)
NE49342	3/4" (19 mm)
NE49344	1" (25.5 mm)
NE49346	1-1/4" (32 mm)

COBB SPINAL ELEVATORS, ALUMINUM HANDLE

All stainless steel working end and anodized aluminum handle. Sharp, 11" (28 cm).




CODE #	DESCRIPTION
NE46413	3/8" (10 mm)
NE46414	1/2" (12 mm)
NE46416	3/4" (19 mm)
NE46418	1" (25 mm)
NE46420	1-1/4" (32 mm)

ENDOSCOPIC ELEVATORS

ENDOSCOPIC RIB ELEVATORS

Sharp, with phenolic handle




Phenolic handles autoclavable to 290°F (143°C)

CODE #	TIP WIDTH	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
NE55640	Angled 45 degrees	9 mm radius	4"	8"	0.45" (11.5 mm)	17" (43 cm)
NE55644	Curved 90 degrees	9 mm radius	4"	8"	0.45" (11.5 mm)	17" (43 cm)

ENDOSCOPIC COBB ELEVATORS

Phenolic handle.




Phenolic handles autoclavable to 290°F (143°C)

CODE #	TIP WIDTH	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
NE55620	1/4" (6 mm)	2.7"	9.3"	0.285" (7.2 mm)	5"	17" (43 cm)
NE55610	1/4" (6 mm)	3-1/2"	10"	0.45" (11.5 mm)	6-1/2"	20" (51 cm)
NE55612	3/8" (10 mm)	3-1/2"	10"	0.45" (11.5 mm)	6-1/2"	20" (51 cm)
NE55614	0.49" (12.5 mm)	3-1/2"	10"	0.45" (11.5 mm)	6-1/2"	20" (51 cm)

SPECIALTY ELEVATORS

DARRACH ELEVATORS


Useful in shoulder surgery. Blunt tips with serrations can be used to elevate or retract soft tissue and muscle from the operative site.



CODE #	DESCRIPTION
NE46902	Small, 1/2" (12 mm) wide at tip, 10-1/4" (26 cm)
NE46903	Medium, 3/4" (19 mm) wide at tip, 10-1/4" (26 cm)
NE46904	Large, 1" (25 mm) wide at tip, 13-3/4" (35 cm)

CHANDLER ELEVATORS

Blunt, slightly curved.



CODE #	TIP WIDTH	LENGTH
NE46250	8 mm	7-1/8" (18 cm)
NE46252	12 mm	8" (20.5 cm)
NE46254	17 mm	9" (23 cm)
NE46256	19 mm	10" (25.5 cm)

HARRINGTON SPINAL ELEVATOR

All stainless steel



CODE #	DESCRIPTION
NE55798	Sharp, 3/4" (19 mm) wide tip, 12" (30.5 cm)

SACHS NERVE ELEVATOR



CODE #	DESCRIPTION
NE57618	3 mm blunt, 7-1/2" (19 cm)