


CURETTES


Curettes are precision instruments used for scraping or debriding bone and soft tissue in orthopedic procedures. Designed with sharp, scoop-shaped tips, they allow for controlled removal of unwanted tissue while preserving surrounding structures. Manufactured from surgical-grade stainless steel, our curettes offer excellent balance, durability, and tactile feedback for precise surgical work.



BRUN CURETTES

BRUN (SPRATT) CURETTES

Oval cup, hollow handle




CODE #	TIP SHAPE	SIZE	LENGTH
NE46511	Straight	#5/0	7" (18 cm)
NE46512	Straight	#4/0	7" (18 cm)
NE46513	Straight	#3/0	7" (18 cm)
NE46515	Straight	#2/0	7" (18 cm)
NE46517	Straight	#0	7" (18 cm)
NE46519	Straight	#1	7" (18 cm)
NE46521	Straight	#2	7" (18 cm)
NE46523	Straight	#3	7" (18 cm)
NE46525	Straight	#4	7" (18 cm)
NE46529	Straight	#5	7" (18 cm)
NE46531	Straight	#6	7" (18 cm)
NE46563	Straight	#3/0	9" (23 cm)
NE46564	Straight	#2/0	9" (23 cm)
NE46565	Straight	#0	9" (23 cm)
NE46566	Straight	#1	9" (23 cm)
NE46567	Straight	#2	9" (23 cm)
NE46568	Straight	#3	9" (23 cm)
NE46569	Straight	#4	9" (23 cm)
NE46570	Straight	#5	9" (23 cm)
NE46571	Straight	#6	9" (23 cm)

CODE #	TIP SHAPE	SIZE	LENGTH
NE49110	Angled 45 degrees	#3/0	9" (23 cm)
NE49111	Angled 45 degrees	#2/0	9" (23 cm)
NE49112	Angled 45 degrees	#0	9" (23 cm)
NE49113	Angled 45 degrees	#1	9" (23 cm)
NE49114	Angled 45 degrees	#2	9" (23 cm)
NE49115	Angled 45 degrees	#3	9" (23 cm)
NE49116	Angled 45 degrees	#4	9" (23 cm)
NE49117	Angled 45 degrees	#5	9" (23 cm)
NE49118	Angled 45 degrees	#6	9" (23 cm)

BRUN (SPRATT) CURETTES


Round up, hollow handle



CODE #	TIP SHAPE	SIZE	LENGTH
NE49120	Straight	#3/0	6-3/4" (17 cm)
NE49121	Straight	#2/0	6-3/4" (17 cm)
NE49122	Straight	#0	6-3/4" (17 cm)
NE49123	Straight	#1	6-3/4" (17 cm)
NE49124	Straight	#2	6-3/4" (17 cm)
NE49125	Straight	#3	6-3/4" (17 cm)
NE49126	Straight	#4	6-3/4" (17 cm)

BRUN CURETTES


8-7/8" (22.5 cm).



CODE #	DESCRIPTION
343500	Size 000
343501	Size 00
343501	Size 0
343503	Size 1
343504	Size 2
343505	Size 3
343506	Size 4

BRUN CURETTES

Solid hexagonal handle, 9" (23 cm).




CODE #	TIP SHAPE	SIZE
NE46551	Straight	#3/0
NE46552	Straight	#2/0
NE46553	Straight	#0
NE46554	Straight	#1
NE46555	Straight	#2
NE46556	Straight	#3
NE46557	Straight	#4
NE46558	Straight	#5
NE46559	Straight	#6

CODE #	TIP SHAPE	SIZE
NE49150	Angled 45 degrees	#3/0
NE49151	Angled 45 degrees	#2/0
NE49152	Angled 45 degrees	#0
NE49153	Angled 45 degrees	#1
NE49154	Angled 45 degrees	#2
NE49155	Angled 45 degrees	#3
NE49156	Angled 45 degrees	#4
NE49157	Angled 45 degrees	#5
NE49158	Angled 45 degrees	#6

COBB CURETTES

COBB CURETTES, ALUMINUM HANDLE

All stainless steel working end and anodized aluminum, knurled handle, 11" (28 cm).




CODE #	TIP SHAPE	SIZE
NE46433	Angled	#3/0
NE46434	Angled	#2/0
NE46435	Angled	#0
NE46436	Angled	#1
NE46437	Angled	#2
NE46438	Angled	#3
NE46439	Angled	#4
NE46440	Angled	#5
NE46441	Angled	#6

CATALOG #	TIP SHAPE	SIZE
NE46422	Straight	#3/0
NE46423	Straight	#2/0
NE46424	Straight	#0
NE46425	Straight	#1
NE46426	Straight	#2
NE46427	Straight	#3
NE46428	Straight	#4
NE46429	Straight	#5
NE46430	Straight	#6

COBB CURETTES, STAINLENE STEEL HANDLE

All stainless steel working end and stainless steel, knurled handle, 11" (28 cm).




CODE #	TIP SHAPE	SIZE
NE46433S	Angled	#3/0
NE46434S	Angled	#2/0
NE46435S	Angled	#0
NE46436S	Angled	#1
NE46437S	Angled	#2
NE46438S	Angled	#3
NE46439S	Angled	#4
NE46440S	Angled	#5
NE46441S	Angled	#6

CODE #	TIP SHAPE	SIZE
NE46422S	Straight	#3/0
NE46423S	Straight	#2/0
NE46424S	Straight	#0
NE46425S	Straight	#1
NE46426S	Straight	#2
NE46427S	Straight	#3
NE46428S	Straight	#4
NE46429S	Straight	#5
NE46430S	Straight	#6

SPECIALTY BONE CURETTES


EPSTEIN CURETTES

Reverse angle, hollow handle, 8" (20.5 cm).



CODE #	DESCRIPTION
NE49580	#000
NE49582	#00
NE49584	#0

SCOVILLE CURETTES



CODE #	DESCRIPTION
NE49574	Straight, 10" (25.5 cm)
NE49576	Down angle, 10" (25.5 cm)
NE49578	Up angle, 10" (25.5 cm)

JANSEN CURETTE



CODE #	DESCRIPTION
NE29632	Double-ended, oval cups, 6" (15 cm)

MARTINI BONE CURETTE



CODE #	DESCRIPTION
NE29620	Double-ended, 4 mm and 5 mm round cups, 5-1/2" (14 cm)

HATFIELD BONE CURETTES

Stainless steel with satin finish. Oval shaped tips, overall length 9" (23 cm). Complete set of nine curettes



CODE #	TIP SHAPE	SIZE	TIP LENGTH
584160	000	2.5 mm	4.0 mm
	Straight	2.8 mm	4.6 mm
	Straight	3.3 mm	5.0 mm
	Straight	3.8 mm	6.0 mm
	Straight	4.8 mm	7.6 mm
	Straight	5.8 mm	9.2 mm
	Straight	7.2 mm	11.6 mm
	Straight	8.7 mm	12.8 mm
	Straight	10.5 mm	15.5 mm

VOLNEANN CURETTES

Double-ended, oval and round cups.



CODE #	DESCRIPTION
NE46560	5" (12.5 cm)
NE46561	8" (20.5 cm)
NE46562	6-1/2" (16.5 cm)


WILLIGER BONE CURETTE



CODE #	DESCRIPTION
NE29630	Double-ended, 3 mm and 4 mm oval cups, 5-1/2" (14 cm)

VOLNEANN CURETTES


Double-ended, oval cups.



CODE #	DESCRIPTION
NE49200	Double-ended, straight, 5" (12.5 cm)
NE49202	Double-ended, straight, 5 1/2" (14 cm)
NE46550	Double-ended, oval cups, 6 1/2" (16.5 cm)
NE49204	Double-ended, straight, 8" (20.5 cm)

RING CURETTES


Hollow handle with fenestrated cup.



CODE #	DESCRIPTION
NE55680	Size 0, small, 10" (25.5 cm)
NE55684	Size 0, small, 11-1/2" (29 cm)
NE55681	Size 2, medium, 10" (25.5 cm)
NE55685	Size 2, medium, 11-1/2" (29 cm)
NE55682	Size 6, large, 10-1/4" (26 cm)
NE55686	Size 6, large, 12" (30.5 cm)

CONE RING CURETTES

Knurled handle and fenestrated cup, 9" (23 cm).



CODE #	TIP SHAPE	SIZE	
NE49450	Straight	#1	8 mm diameter
NE49452	Straight	#2	6 mm diameter
NE49454	Straight	#3	3 mm diameter
NE49456	Angled 35 degrees	#1	8 mm diameter
NE49458	Angled 35 degrees	#2	6 mm diameter
NE49460	Angled 35 degrees	#3	3 mm diameter

RAY PITUITARY CURETTES

7-1/2" (19 cm).




CODE #	TIP SHAPE	SIZE
NE49403	Vertical angle 45 degrees	3 mm
NE49404	Vertical angle 45 degrees	4 mm
NE49405	Vertical angle 45 degrees	5 mm
NE49413	Vertical angle 90 degrees	3 mm
NE49414	Vertical angle 90 degrees	4 mm
NE49415	Vertical angle 90 degrees	5 mm
NE49423	Horizontal angle 45 degrees	3 mm
NE49424	Horizontal angle 45 degrees	4 mm
NE49425	Horizontal angle 45 degrees	5 mm
NE49433	Horizontal angle 90 degrees	3 mm
NE49434	Horizontal angle 90 degrees	4 mm
NE49435	Horizontal angle 90 degrees	5 mm

ENDOSCOPIC CURETTES


ENDO BRUN CURETTES

Endoscopic with phenolic handle



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
NE55532	Angled 40°	2 x 4 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
NE55523	Straight	3 x 6 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
NE55533	Angled 40°	3 x 6 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
NE55543	Angled 40°	3 x 6 mm cup	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)
NE55556	Angled 30°	6 x 9 mm cup	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)


ENDO EPSTEIN CURETTES



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
NE55593	Reverseangled 40°	3 x 6 mm	2.7"	9.3"	0.285"	5"	17" (43 cm)
NE55586	Reverseangled 30°	6 x 9 mm	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)

ENDO SCOVILLE CURETTE


Endoscopic with phenolic handle.



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
NE55510	Straight	9 x 15 mm	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)

ENDO CONE RING CURETTES

Endoscopic with phenolic handle



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
NE55573	Angled 40°	3 mm	2.7"	9.3"	0.285"	5"	17" (43 cm)
NE55566	Angled 30°	6 mm	10"	10"	0.45"	6-1/2"	20" (51 cm)

SMARTRACK CURETTE TRAYS

BRUN SMARTRACK FOR 7" CURETTES

Aluminum and stainless steel.



CODE #	DESCRIPTION
NE98160	Holds 12 7" NEedic® Brun Cures

COBB SMARTRACK FOR 11" INSTRUMENTS

Aluminum and stainless steel.



CODE #	DESCRIPTION
NE98170	Holds nine 11" NEedic® Cobb Cures, Elevators or Gouges

HIBBS AND BRUN SMARTRACK FOR 9" INSTRUMENTS

Aluminum and stainless steel.



CODE #	DESCRIPTION
NE98180	Holds 12 9" NEedic® Brun Cures, or Hibbs Chisels, Gouges or Osteotomes