



Our orthopedic screwdrivers are designed for precise control and secure handling of all standard implant screws. With ergonomic handles and torque-optimized shafts, they allow for accurate insertion and removal without stripping or slippage. Available in a variety of tip designs and sizes, each screwdriver ensures compatibility with a wide range of orthopedic systems.

SCREWDRIVERS

HEX HEAD HEXAGONAL SCREWDRIVERS



CODE #	DESCRIPTION
NE47952	2.5 mm screwdriver with sleeve, 7-7/8" (20 cm)
NE48368	2.5 mm screwdriver with notch, 8-1/4" (21 cm)
NE48349	3.5 mm screwdriver with notch, 9-5/8" (24.5 cm)
NE47955	Screw holding sleeve only for NE47952, 2" (5 cm)
NE48365	Screw holding sleeve only for NE48368, 3" (8 cm)
NE47983	Screw holding sleeve only for NE48349, 5" (12.5 cm)
NE47956	Screwdriver only for NE47952

HEXAGONAL SCREWDRIVER SHAFT (PELVIC)




CODE #	DESCRIPTION
NE47951	2.5 mm with notch, quick coupling end, 6-1/2" (16.5 cm)
NE47953	2.5 mm without notch, quick coupling end, 4" (10 cm)
NE47961	3.5 mm with notch, quick coupling end, 6-1/2" (16.5 cm)
NE48365	Sleeve only for NE47951, 3" (8 cm)
NE47983	Sleeve only for NE47961, 5" (12.5 cm)

SCREWDRIVER SLEEVES



CODE #	DESCRIPTION
NE48365	Sleeve only for NE47951, NE48348, NE48364, NE48368, 3" (8 cm)
NE47983	Sleeve only for NE48349, NE48353 and NE47961, 5" (12.5 cm)

LARGE HEXAGONAL SCREWDRIVER (PELVIC)




CODE #	DESCRIPTION
NE48353	3.5 mm screwdriver with notch, 11-3/4" (30 cm)
NE48364	Screwdriver only, 2.5 mm with notch, 11" (27 cm)
NE47983	Sleeve only for NE48353, 5" (12.5 cm)

LARGE HEXAGONAL SCREWDRIVER



CODE #	DESCRIPTION
NE48252 3	3.5 mm, quick coupling end, 3-7/8" (10 cm)

LARGE HEX SCREWDRIVER WITH UNIVERSAL JOINT



CODE #	DESCRIPTION
NE48396	3.5 mm tip, 11-5/8" (29.5 cm)

HEXAGONAL SCREWDRIVER SHAFT



CODE #	DESCRIPTION
NE48348	2.5 mm with notch, quick coupling end, 5-1/2" (14 cm)
NE48365	Sleeve only

LARGE HEXAGONAL SCREWDRIVER



CODE #	DESCRIPTION
NE47954	3.5 mm, quick coupling end, 3-7/8" (10 cm)

SCREWDRIVER




CODE #	DESCRIPTION
NE48670	1.5 mm hexagonal, 7-7/16" (200 mm)
NE48350	2.5 mm hexagonal, 9-27/32" (250 mm)
NE48352	3.5 mm hexagonal, 9-27/32" (250 mm)
NE48354	Single slot, 9-27/32" (250 mm)
NE48356	Cruciform, 9-27/32" (250 mm)
NE48358	Phillips, 9-27/32" (250 mm)

MINI HEXAGONAL SCREWDRIVER



CODE #	DESCRIPTION
NE47950	1.5 mm, mini quick coupling end, sleeve included, 2-1/4" (5.5 cm)
NE47950B	Shaft only, 1.5 mm with notch, mini quick coupling end, 2" (5 cm)
NE47959	Screw holding sleeve only

SINGLE SLOT SCREWDRIVER



CODE #	DESCRIPTION
NE48355	Phenolic handle, 6-3/4" (17 cm)

CRUCIFORM HEAD CRUCIFORM SCREWDRIVER



CODE #	DESCRIPTION
NE47957	1.5 / 2.0 mm screwdriver with sleeve, 7-1/2" (19 cm)
NE47959	Screw holding sleeve only for NE47957
NE47958	Screwdriver only for NE47957

MINI CRUCIFORM SCREWDRIVER BIT



CODE #	DESCRIPTION
NE45004	1.5 / 2 mm, mini quick coupling end, 2-1/2" (6.5 cm)

T-HANDLE LARGE HEX SCREWDRIVER W/ T-HANDLE

CODE #	DESCRIPTION
NE48395	3.5 mm, quick coupling end, 6-7/8" (17.5 cm)

DEPTH GAUGE SCREW DEPTH GAUGES

Used to determine screw length required during internal fracture fixation (mm graduated). Snaps apart for easy cleaning. IMPORTANT NOTE: When not in use, retract measuring end inside gauge for probe protection



CODE #	DESCRIPTION
NE47934	Mini, up to 30 mm, 4-3/4" (12 cm)
NE47936	Small, up to 50 mm, 6-1/4" (16 cm)
NE47940	Large, up to 100 mm, 8-3/4" (22 cm)
NE47935	Extra-long, up to 120 mm, 11" (28 cm)

INTERCHANGEABLE SCREWDRIVER SET

Screwdriver Sets



CODE #	DESCRIPTION
NE48473	Screwdriver handle. Allows K-Wire and pins to be inserted with Jacobs Universal Chuck NE48478. Cannulation 1/4" (6.3 mm), includes key, 5-1/2" (14 cm)
NE48474	Single slot, 5-1/2" (14 cm)
NE48475	Cruciform, 5-1/2" (14 cm)
NE48476	Phillips, 5-1/2" (14 cm)
NE48477	Small Phillips, 5-1/2" (14 cm)
NE48478	Jacobs Universal Chuck. Allows K-Wire and pins to be inserted with Jacobs Chuck. Cannulation 1/4" (6.3 mm), includes key, 3" (7.5 cm)
NE48478P	Key for NE48478 Jacobs Universal Chuck
NE48479	3.5 mm diameter, 5-1/2" (14 cm)
NE48480	2.5 mm diameter, 5-1/2" (14 cm)

SCREW HOLDING FORCEPS

SELF-RETAINING SCREW HOLDING FORCEPS



CODE #	DESCRIPTION
NE47930	3-1/2" (9 cm)

SELF-RETAINING SCREW HOLDING FORCEPS



CODE #	DESCRIPTION
NE45208	For 1.5 mm to 2.7 mm screws
NE45210	For 3.5 mm to 6.5 mm screws